

Table SI. Sequence of primers for limbic system-associated membrane protein break-point verification.

Case	Forward (5'-3')	Reverse (5'-3')	Amplicon size, bp
NBP2T	AATTGGAATGCTTAAAGGA CATTC	TAGTGTACTACTAGGCAGT CAG	438
NBP5T (1)	AGGGACTTAAGTTATAAAA GAACAG	TGTCAGCTATATGTAACAT GG	653
NBP5T (2)	GAAACCCCGTCTCTACTAAA TG	TCCTTCTTAACCTCCCTAA CC	237
NBP5T (3)	GACAGTCTTGATTTCTCAC C	CACCCTCTCATTTTATTCA AAAGCC	276
NBP10 T	TGTCTAGGGCATGAGCGGA G	ATTCATACACTGGAGGTCA CAGG	300
NBP14 T	AAGATGCTCATTCCCCTCC	TCCACTTTTGTGTGCTAAC C	247
NBP26 TR	GGGACATATTGAAAGTGCT GAG	CAGACACATGGAACTACTG C	285
NBP28 TR	GCGAACATTTACATTTCACT TCCC	GCTATAAAAGCAAGAGCA AACC	307