

Figure S1. Sequencing primers used for BRAF V600E mutation testing. (A) Forward and (B) reverse BRAF primers. (C) Forward and (D) reverse M13 primers. F, forward; R, reverse.

A	BRAF-F	5'- TGTAAAACGACGGCCAGTGACCTTCAATGACTTTCTAGTAACTCAGC -3'
B	BRAF-R	5'- CAGGAAACAGCTATGACCACCTAAACTCTTCATAATGCTTGCTCT -3'
C	M13-F	5'- TGTAAAACGACGGCCAGT -3'
D	M13-R	5'- CAGGAAACAGCTATGACC -3'

Figure S2. Sequencing traces for Sanger sequencing.  
(A) Positive and (B) negative for the BRAF V600E mutation.

